Bodem – Bepaling van ammoniumstikstof

1 PRINCIPE

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3.2 REAGENTIA

- 3.2.1 Kaliumchloride-oplossing, 1 mol/l: 74.6 g/l KCl in water

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- ISO/TS 14256-1:2003 Soil quality - Determination of nitrate, nitrite and ammonium in field-moist soils by extraction with potassium chloride solution - Part 1: Manual method
- ISO 14256-2:2005 Soil quality - Determination of nitrate, nitrite and ammonium in field-moist soils by extraction with potassium chloride solution - Part 2: Automated method with segmented flow analysis
- NBN EN ISO 11732:2005 Water quality — Determination of ammonium nitrogen — Method by flow analysis (CFA and FIA) and spectrometric detection (ISO 11732:20005)
- NBN EN ISO 15923-1:2024 Water quality - Determination of selected parameters by discrete analysis systems - Part 1: Ammonium, nitrate, nitrite, chloride, orthophosphate, sulfate and silicate with photometric detection (ISO 15923-1:2013)**

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De rapportagegrens is $\leq$ **XX** kg $\text{NH}_4\text{-N/ha}$ per bodemlaag.

8 REFERENTIES

- a. ISO/TS 14256-1:2003 Soil quality - Determination of nitrate, nitrite and ammonium in field-moist soils by extraction with potassium chloride solution - Part 1: Manual method

- b. ISO 14256-2:2005 Soil quality - Determination of nitrate, nitrite and ammonium in field-moist soils by extraction with potassium chloride solution - Part 2: Automated method with segmented flow analysis
- c. NBN EN ISO 11732:2005 Water quality — Determination of ammonium nitrogen — Method by flow analysis (CFA and FIA) and spectrometric detection (ISO 11732)
- d. NBN EN ISO 15923-1:2024 Water quality - Determination of selected parameters by discrete analysis systems - Part 1: Ammonium, nitrate, nitrite, chloride, orthophosphate, sulfate and silicate with photometric detection (ISO 15923-1:2013)
- e. NEN 6604:2007 Water - Bepaling van het gehalte aan ammonium, nitraat, nitriet, chloride, ortho-fosfaat, sulfaat en silicaat met een discreet analysesysteem en spectrofotometrische detectie.